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LaboPol-11-21-41-51-6
Instruction Manual

Technical Data

Subject		Specifications				
		LaboPol-1	LaboPol-2	LaboPol-4	LaboPol-5	LaboPol-6
Disc	Diameter	200-230mm / 8-9"				
	Speed	250 rpm	250/500 rpm	0-150 rpm	50-500 rpm	100-1200 rpm
	Rotational direction	Counter-clockwise				
	Motor	180 W	90/180 W	120 W	370 W	550 W
	Torque at disc Continuous Max.	6.9 Nm 12.6 Nm	3.5 Nm 8.6 Nm	5.7 Nm 11.0 Nm	5.2 Nm 10.4 Nm	5.0 Nm 10.0 Nm
Directives and Standards	General 98/37/EC & 73/23/EEC					
	EMC 89/336/EEC & 92/31/EEC					
Noise level	at idle running at a distance of 1.0 m/39.4" from the machine	Approx. 58 dB(A)	Approx. 58 dB(A)	Approx. 53 dB(A)	Approx. 65 dB(A)	Approx. 65 dB(A)
Surrounding temperature		5-40°C / 41-104°F				
Humidity	Non condensing	0-95%RH				
Supply	Power	1 (N+L+PE) phase				
	Power consumption	270 W	245 W	230 W	580 W	750 W
	Voltage / frequency	Max. continuous load				
	100V / 50 Hz	3.0 A	3.0 A			
	100V / 60 Hz	2.8 A	2.9 A			
	100-120V / 60 Hz	2.7 A	2.6 A			
	220-240V / 50 Hz	1.3 A	1.3 A			
	220-240V / 60 Hz	1.2 A	1.3 A			
	100-120V / 50-60 Hz			2.5 A	5 A	7.1 A
	220-240V / 50-60 Hz			1.25 A	2.4 A	3.6 A
Pressure for water tap	1-10 bar / 14.5-145 psi					
Water inlet	1/2" or 3/4"					
Water outlet	ø32mm / 1 1/4"					
Dimensions and weight	Width	410 mm / 16.1"				
	Depth	670 mm / 26.4"				
	Height	315 mm / 12.4"				
	Weight	22 kg / 48.5 lbs	28 kg / 61.7 lbs	22.5 kg / 49.6 lbs	22.5 kg / 49.6 lbs	27.5 kg / 60.6 lbs

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

*Exotom / Exotom-M
Spare Parts and Diagrams*


Fig. 13


Fig. 14

*Exotom / Exotom-M
Spare Parts and Diagrams*

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
Hot Mounting Resins

Resin	Application	Specific property	Material
ConduFast	Electropolishing	Electrically conductive Very low shrinkage	Acrylic resin with iron powder filler, Thermoplastic, Gray
SpeciFast	Glass clear mounts Porous samples Surface electrical insulator for ConduFast	Transparent Medium shrinkage	Acrylic Thermoplastic
PolyFast	Edge retention	Very low shrinkage Medium removal rate	Phenolic resin with carbon filler, Thermosetting, Black
	Examination in scanning electron microscope	Electrically conductive Low emission in the vacuum chamber	
IsoFast	Edge retention Planeness	Low shrinkage Low removal rate	Diallylphthalate with glass fibre filler, Thermosetting, Green
DuroFast	Edge retention Planeness	Very low shrinkage Good adhesion Very low removal rate	Epoxy with mineral filler, Thermosetting, Black
MultiFast	Routine examination Backup resin	Low shrinkage Medium removal rate	Phenolic resin with wood flour filler, Thermosetting, Black
MultiFast Green	Routine examination Colour marking	Low shrinkage Medium removal rate	Phenolic resin with wood flour filler, Thermosetting, Green
MultiFast Brown	Routine examination Colour marking	Low shrinkage Medium removal rate	Phenolic resin with wood flour filler, Thermosetting, Brown
Pre-Mount	Serial mounting of uncomplicated shapes	Easy to handle Low shrinkage Medium removal rate	Phenolic resin with wood flour filler, In tablet form, Thermosetting, Black


Curing times	Dia. 25 mm (1")	Dia. 30 mm (1¼")	Dia. 40 mm (1½")	Dia. 50 mm
Heating / Cooling time	min.	min.	min.	min.
ConduFast	4 / 3	5 / 3	6 / 4	7 / 4
SpeciFast	5 / 8	6 / 8	7 / 10	8 / 12
PolyFast	5 / 3	6 / 3	7 / 4	8 / 4
IsoFast	5 / 3	6 / 3	7 / 4	9 / 4
DuroFast	5 / 3	6 / 3	7 / 4	8 / 4
MultiFast	5 / 3	6 / 3	7 / 4	9 / 4
MultiFast Green*	7 / 3	8 / 3	10 / 4	12 / 4
MultiFast Brown	5 / 3	6 / 3	7 / 4	9 / 4
Pre-Mount	5 / 3	6 / 3	7 / 4	9 / 4

The heating times refer to a heating temperature of 180 °C. The heating and cooling times may be shorter for samples with large surface areas and high heat conductivity.
* MultiFast Green has longer heating times due to a lower heating temperature of 150 °C

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Cold Mounting Materials


Name	Application	Features	Curing time	Mounting cup
Specifix-20	Epoxy	Vacuum impregnation Porous samples Mineralogical samples	8 h.	Multiform
Specifix-40	Epoxy	Vacuum impregnation Porous samples Mineralogical samples	3.5 h. at 50 °C	Multiform
Epofix	Epoxy	Vacuum impregnation Porous samples Mineralogical samples	8 h.	Multiform
Caldofix	Epoxy	Vacuum impregnation Porous samples Mineralogical samples	1.5 - 2 h. at 70-80 °C	Multiform
MetaFix-20	Acrylic	Serial mounting Irregularly shaped samples	15 min.	Multiform Flexiform Seriform
Citofix	Acrylic	Serial mounting Irregularly shaped samples	7 - 10 min.	Multiform Flexiform Seriform
Durofix-2	Acrylic Mineral filler	Serial mounting Edge retention Irregularly shaped samples	15 min.	Multiform Flexiform Seriform
Acryfix	Acrylic	For routine work	9 - 10 min.	Multiform Flexiform Seriform
Triofix-2	Polyester / Acrylic / Mineral filler	Edge retention Planeness	15 - 18 min.	Multiform Flexiform
Serifix	Polyester	Serial mounting	45 min.	Multiform Flexiform Seriform

Carl Zeiss

INTRODUCTION
Overall view

Photomicrography using
35 mm SLR cameras

Overall View of 35 mm SLR Camera on 63 mm Objective and Adapter for Microscope Camera (Image Position 44 mm) with 10x Camera Eyepiece on 44 mm Camera Port – i.e. on the Axioskop with Phototube in this Example

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
Photomicrography using
35 mm SLR cameras

INTRODUCTION
Overall view

Carl Zeiss

Overall View of 35 mm SLR Camera on 2.5x / T2 Adapter on Front Camera Port of the Axiovert 25 C

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